

Archiving cultures at $-80\text{ }^{\circ}\text{C}$

Adapted from the lab compendium Snowflake Yeast Protocols — Ratcliff Lab / MuLTEE (V2).

Archiving for storage in $-80\text{ }^{\circ}\text{C}$ freezer

1. In a 1.5 mL tube, combine **750 μL culture** with **750 μL 70% glycerol** (sterile), for a final glycerol concentration of 35%.
2. Vortex at maximum speed ($\sim 3,000$ rpm) briefly.
3. Label the tube (population, treatment, timepoint/transfer number, date) and transfer immediately to $-80\text{ }^{\circ}\text{C}$.