

Daily growth and settling selection

Adapted from the lab compendium Snowflake Yeast Protocols — Ratcliff Lab / MuLTEE (V2).

Growth

1. Inoculate from the previous day's settled fraction (see [[> Settling selection](#)]) or from a fresh colony/stock into 10 mL fresh growth medium (YPD, per liter: 1 L water, 20 g dextrose, 20 g peptone, 10 g yeast extract, autoclaved).
2. Grow with shaking for approximately 24 h, at **30 °C, 250 rpm**.

Settling selection

1. Vortex the culture very briefly to resuspend clusters. Do not vortex at over **~2,300 rpm**.
2. Transfer 1 mL of culture into a new tube to settle. Use a wide-bore pipette tip if working with macroscopic clusters.
3. Settle undisturbed:
 - a. Macroscopic populations: 30 seconds.
 - b. Microscopic populations: 3 minutes.
4. Remove the top 900 μL by eye without disturbing the settled pellet (pipette from above the 0.1 mL mark on the tube).
5. Transfer 50 μL from the bottom of the remaining $\sim 100 \mu\text{L}$ into fresh medium to start the next day's culture (using a wide-bore pipette tip if working with macroscopic clusters).