

Reviving strains from frozen stock

Adapted from the lab compendium Snowflake Yeast Protocols — Ratcliff Lab / MuLTEE (V2).

Reviving strains from frozen stock

1. Using a P10 pipette tip, scrape a small amount of ice directly from the surface of the frozen stock.
2. Inoculate the scraped material into **10 mL fresh YPD**.
3. Grow out following the standard growth protocol ([[Daily growth and settling selection > Growth]]) before beginning any re-acclimation or further propagation.
4. Go through 3-5 rounds of settling selection ([[Daily growth and settling selection > Settling selection]]) to encourage growth back to full size. Taking 15 μ L transfers instead of 50 μ L transfers from settled culture sometimes accelerates this process.